

25/5/22

Capillary Electrophoresis.

Slab/plate electrophoresis:

- ① It is slow, labour intensive, difficult to automate and does not yield definitive quantitative results.

- ★ ① Capillary electrophoresis : is an instrumental version and a substitute for slab electrophoresis.
- ② It displays high-speed performance, high resolution separation, requires small sample size i.e., 0.1 to 10 nanometres.
- ③ Quantitative detectors can be used instead of staining.

Instrumentation.

Requirement:


- ① Buffer filled fused silica capillary, 10 to 100 micrometres for internal diameter and 40-60 cm long.
- ② Two platinum electrodes.
- ③ High Voltage supply.
- ④ Sample injector.

⑤ Sample inlet and outlet vial

⑥ Detector or

⑦ Buffer Solution.

capillary



* sample injection

① hydrodynamic injection

a) by applying pressure.

b) by applying vacuum.

c) by gravitation.

② Electro kinetic injection: ie; by using electric supply.

Detectors:

- ① UV-vis absorption spectrometric detectors.
- ② Other detectors include fluorescence, amperometric, electrochemiluminescence detectors.
- ③ Mass spectrometric.

~~A fiber capillary result~~

Electrophogram

It is similar to chromatogram plot of time from injection on x axis v/s the detector signal on y axis.

Working:

- ① Components include sample ^{vial}, source and destination vial, capillary, electrodes, data output and handling devices.
- ② source vial, destination vial and capillary are filled with an electrolyte (buffer solution).
- ③ To introduce the sample, capillary inlet is placed into the sample vial.
- ④ Sample is introduced in capillary by a

Capillary action, pressure, or Electrokinetically and then capillary is returned to source vial.

- ③ Migration of Analites is initiated by electric field applied between the source and destination vials and the is supplied through electrodes by high voltage supply.
- ④ Often in CE all the ions be it negative or positive are pulled through the capillary in the same direction because of electro osmotic flow (EOF).
- ⑤ Analites separated as they migrate due to electrophoretic mobility and are detected using detectors near the outlet end of capillary.
- ⑥ Output send to data output handling device such as integrator or computer.
- ⑦ Data displayed as electropherogram.
- ⑧ Separated compounds appear as peaks with different retention time in electropherogram.
- ⑨ Retention time is the time taken for a solute to pass through the column

Applications:

- 1) Illicit drug analysis
- 2) Analysis of GSR and explosive constituents.
- 3) To investigate in composition.
- 4) Analysis of small molecule pharmaceuticals and in their preparations.
- 5) Presence of pharmaceutical products in standard solutions and biological fluids.